



Mobile hand wash station



HMW230 The MediWash

Matches World Health Organization Guidelines
"Clean hands protect against infection"

The automatic sensor starts a pre-programmed hand wash cycle when hands are placed within range.



Handwash stations must "...be used for that purpose only" - National infection prevention and control manual for England

The MediWash can be plugged in to any convenient 13 amp socket

Designed for medical and healthcare settings

 Order online at tealwash.com

The MediWash mobile handwash station is entirely self contained and needs no mains water service or drainage connections.

Pre-programmed hand wash cycle

The MediWash wash / rinse cycle is designed to match the World Health Organization Guidelines "Clean hands protect against infection".

To start the timed wash programme:-

1. Place the hands by the automatic sensor and it will flash.
2. Water runs for 4 seconds to enable the wetting of hands.
3. The water stops for 20 seconds to allow for the application of liquid soap and performance of the 9 step hand wash routine.
4. Appropriate hand washing techniques are displayed on the screen.
5. Hot water runs for 10 seconds to enable the soap to be rinsed off the hands.
6. The total wash time is approximately 35 seconds.

MediWash features:-

- Hands free operation - hot water hand washing without the need for external plumbing or drainage.
- Plug in to any convenient 13 amp socket.
- Nothing to reset or adjust, with a simple daily hygiene routine.
- The waste water is contained within a heavy duty polyethylene bag - when the fresh is used up the waste bag will be full.
- Up to 40 wet/soap/rinse cycles per filling of each of the 9 litre TEALTainers.
- *Please note:-* the MediWash has nothing to adjust. Both the water temperature (approx. 38 - 40°C) and the wet/soap/rinse cycles are factory set.

Specification

- Overall height 1370 mm
- Height to bowl 880 mm
- Width 620 mm
- Depth 590 mm
- Weight (dry) 30.8 kg (without dispensers)
- TEALTainer weight 1.5 kg (empty), 11.5 kg (full)
- Power supply 230 V mains

Construction

- The body is fully moulded in blue fleck high impact polyethylene.
- The bowl is white acrylic capped ABS.
- Comes complete with 13 amp moulded plug.
- TEALTainers securely retain waste water in a separate heavy duty LLDPE polythene bags.

Included with the MediWash mobile handwash station

- 2 TEALTainer water containers.
- The backboard is designed to accept the customers' own choice of towel and soap dispenser.



Handwash stations must "...be used for that purpose only" - National infection prevention and control manual for England

 Order online at tealwash.com



Teal Patents Ltd

Unit 2
Chelmsley Wood Industrial Estate
Waterloo Avenue
Chelmsley Wood
Birmingham, B37 6QQ, UK
Tel: +44 (0) 121 770 0593
Email: enquiries@tealwash.com

Order online at www.tealwash.com